FLUSHTECHG



- Pulsating Flush Technology
- Flushing Sensor with Alarm
 - Manual or Auto Flush Settings
- High Capacity Flushing
 - Reporting Capabilities
- No Syringing Required
 - System Decontamination Cycle



Automated Endoscope Flushing No Manual Syringing

The FlushTech GI is a powerful, easy-to-operate endoscope flushing aid that uses single and dual flushing modes to move chemistry and water through all of the channels of an endoscope. Pulsating flush technology works to help remove bio-burden from the endoscope

before high-level disinfection in an automated endoscope reprocessor. It includes built-in programmable settings, allowing users to set flushing parameters that will meet specific conditions for achieving optimum cleaning results and high capacity processing.



Product Specifications

Specification	Description
Dimensions	7" H x 8" W x 7" D (17.8 cm x 20.3 cm x 17.8 cm)
Case Rating	IP-53
Flush Pump Flow Rate	.75 GPM (2.84 LPM)
Flush Pump Suction	Maximum suction height - 6ft. (1.82m)
Power Supply/Voltage	Wall Mount Type, In: 110-240 VAC 1A 50-60 Hz, Out: 24 VDC 1.5A
Chemical Compatibility	Flush Pump – Industry standard high level disinfectant, enzymatic and other detergents in diluted form. Do not use Isopropyl Alcohol with the flush pump.
Flush Tubes/Hook Up's Compatibility	Industry standard enzymatic and other detergents in diluted form. Do not use Isopropyl Alcohol with the flush pump.
Unit Weight	9.4 lbs (4.3 Kg)
Approvals	Pending
	The state of the s


Ordering Information

Part Number	Description
7260030	Knight FlushTech GI with power supply adapters for North America,
	North Asia, Europe, UK, China, Australia

Repor	t Sample	е		
Serial #	Date	Time	Event Time	Event Type
0000000	11/15/2013	03:16:26PM	0075	Decontamination
3298567	11/15/2001	03:16:26PM	0010	Single Channel
3298567	11/15/2013	03:16:26PM	0010	Single Channel
3298567	11/15/2013	03:16:26PM	0010	Single Channel
3875649	11/15/2013	03:16:26PM	0010	Single Channel
3875649	11/15/2013	03:16:26PM	0010	Single Channel
3875649	11/15/2013	03:16:26PM	0010	Single Channel
2783645	11/15/2013	03:16:26PM	0010	Single Channel
2783645	11/15/2013	03:16:26PM	0010	Single Channel
2783645	11/15/2013	03:16:26PM	0010	Single Channel
2836243	11/15/2013	03:16:26PM	0010	Single Channel
2836243	11/15/2013	03:16:26PM	0010	Single Channel
2836243	11/15/2013	03:16:26PM	0010	Single Channel
0000000	11/15/2013	03:21:13PM	0075	Decontamination
3298567	11/16/2013	03:21:49PM	0010	Dual Channel
3298567	11/16/2013	03:21:49PM	0010	Dual Channel
3298567	11/16/2013	03:21:49PM	0010	Dual Channel

Accessories

The FlushTech GI comes with a complete set of endoscope irrigation flushing tube components for the most common flexible endoscopes.





A Unit of IDEX Corporation 20531 Crescent Bay Drive Lake Forest, CA 92630 USA

Toll Free: (800) 854-3764 Tel: (949) 595-4800 Fax: (949) 595-4801 www.knighthc.com **Knight Canada**Tel (905) 542.2333
Fax (905) 542.1536

Knight North Asia Tel +82 (2) 3481.6683 Fax +82 (2) 3482.5742 **Knight Europe** Tel +44 (1323) 514.855 Fax +44 (1323) 514.828

Knight South Asia Tel +65 9170.0984 Fax +65 6489.6723

Knight Australia Tel +61 (02) 9352.1801 Fax +61 (02) 9352.1899